

Hons practical MC 4

GROUPS:

1A

- Anina Rademeyer
- Craig October
- Linique Mc Callum

2B

- Molatelo Peter
- Karli Bothma
- Hylton Gibson

3C

- Carmen de Villiers
- Bianca Nilson
- Ankia Visser

4D

- Melissa De Lilly
- Lee-Maine Spies
- Jessica Reid

5A

- Elanda Relling
- Tyron Nel
- Christopher Borrageiro

6B

- Andy van der Berg
- Claire du Plessis
- Sandisiwe Matyesini

7C


- Lieke Dale
- Amy Hare
- Sarah Kellow-Web

8D

- Nicolaas Grobler
- Aqilah Benjamin
- Aneliswe Ngcobo

ABOUT THE PRAC:

(1) Characterise the enzyme:



(2) The experiments:

GROUP A

- Reverse direction
- Varying [NAD⁺]
- Same [Lac] +[LDH] + [Hydrazine]

GROUP B

- Reverse direction
- Varying [Lac]
- Same [NAD⁺] +[LDH] + [Hydrazine]

GROUP C

- Forward direction
- Varying [NADH]
- Same [Pyr] +[LDH]

GROUP D

- Forward direction
- Varying [Pyr]
- Same [NADH] +[LDH]

GROUP A

- Product Inhibition: NADH
- Varying [NADH]
- Same [Lac] +[NAD⁺] + [LDH] + [Hydrazine]

GROUP B

- Without Hydrazine
- Varying [Lac]
- Same [NAD⁺] +[LDH]

GROUP C

- Product Inhibition: NAD⁺
- Varying [NAD⁺]
- Same [Pyr] +[NADH]+[LDH]

GROUP D

- Product Inhibition: Lac
- Varying [Lac]
- Same [NADH] + [Pyr]+[LDH]

(3) An enzyme assay:

- Pipette metabolites into the cuvettes, the reaction is initiated by the addition of the enzyme.
- Once the enzyme is added measure the absorbance readings for the reaction over a period of 3mins.
- The types of data we expect to see- Forward reaction: decreasing Abs; Reverse reaction: increasing Abs

(4) What needs to be done:

Create a pipetting schema so that we investigate 10 different concentrations around the Km values. Take the following into consideration:

a- Given :

Km	Stock concentrations
Km(NADH) = 0.05 mM	[NADH] = 1.5 mM
Km(NAD ⁺) = 0.4 mM	[NAD] = 5 mM
Km(Pyr) = 0.5 mM	[Pyr] = 10 mM
Km(Lac) = 10 mM	[Lac] = 1 M

HYDRAZINE is given at a 10X dilution of bottle concentration

b- Volumes:

- The total cuvette volume is 1 ml and of that you use 10 µl of the enzyme (LDH) and if you are doing **reverse direction experiments (only)** 35 µl is the Hydrazine.
- The other volumes are of the substrate and co-factor which you calculated in the pipetting schema.
- Lastly, the remaining volume to make up this 1 ml volume is done by adding buffer to the cocktail mixture.

c- How to determine the 10 concentrations:

- The range you want to investigate is between: [Km/10; Km] and (Km; 10Km] and the fixed substrates are usually at saturating condition i.e. 10xKm. Use Excel to create an incremental pipetting schema.

(5) The Program:

14/03: Wednesday	15/03 Thursday:	16/03 Friday:
Morning: Practical theory and pipetting schema	Morning: 09H30 Research lecture	Morning: Group 1-4 (second set of experiments)
Afternoon: Group 1-4 (first set of experiments)	Mid Morning: Group 5-8 (first set of experiments)	Afternoon: Group 5-8 (second set of experiments)

19/03 Monday: Redo experiments Afternoon: Data analysis	20/03 Tuesday: Morning: Data analysis Afternoon: Seminar	21/03 Wednesday: Human rights day
22/03 Thursday: Lecture (9:30) Data analysis (10:30)	23/03 Friday: Data analysis	26/03 Monday: Hand in report (16:00)